

**This product may be purchased from Connevens Limited secure online store
at www.DeafEquipment.co.uk**



DeafEquipment.co.uk

Solutions to improve the quality of life



The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

Part No. 10298869 / JN4468 / Sept. 07
Printed in the UK

Siemens Hearing Instruments Ltd.
Alexandra House
Newton Road
Manor Royal
Crawley RH10 9TT

Tel: 01293 423700
Fax: 01293 403080
email: siemenshearing.shi.uk@siemens.com

www.siemens.co.uk/hearing



The HearCheck Screener

The easy, reliable way to triage patients with hearing loss

www.siemens.co.uk/hearing

SIEMENS



Fast, accurate, simple and safe identification of hearing loss: the HearCheck Screener

Studies have shown that hearing loss frequently goes undetected – and untreated

Fast, accurate, simple and safe

For both patients and testers alike, establishing hearing loss with the technologically advanced HearCheck Screener is fast, accurate, simple and hygienic. The test can be performed in a mere 30 seconds per ear and the results are available immediately.

What is more, with its one-button operation, the HearCheck Screener is extremely easy to use, requiring no special skills.

It features two test frequencies 1kHz, 3kHz with different levels in the range of 20 to 75 dBHL. The HearCheck Screener results can be used to streamline referrals for full audiometric assessment and one of today's advanced range of digital hearing systems such as open / instant fit or custom earmould models.

Enhanced benefits

Patients, audiologists and dispensers can all benefit from increased referral efficiency at the primary care level.

Testing in four easy steps

Testing with the HearCheck Screener couldn't be easier, involving just four simple steps:

- 1) Gently hold the device to the patient's ear (make sure the ear cover completely surrounds the skin surface around the ear)
- 2) Instruct patient to indicate when they hear a tone
- 3) Press the Start button once to initiate test sequence (actual testing begins after completion of a short functional test, indicated by three flashing lights)
- 4) Record the results



The HearCheck advantage

- Improves hearing assessment efficiency
- Supports fitting decision-making
- Easy LED colour-coded assessment
- Easy to use with quick identification of degree of hearing loss.
- CE approved medical device class 2a

Technical details

- Two test frequencies (1kHz, 3kHz)
- Available test levels: 20, 35, 55, 75 dBHL
- Accuracy: 3 dB
- Easy-to-understand LED colour-coded testing assessment (green, amber, red)

HearCheck Screener Package

- One device per box
- Starter pack of ear cups (25 pcs.)
- 2 AAA size batteries
- User manual