

# Personal stereo direct input lead

The blue plugged personal stereo direct input lead allows direct input hearing aids to connect directly to a 3.5mm headphone socket.

Suitable for use with iPods, portable CD players, computers, keyboards etc. This lead acts as a direct replacement for a pair of headphones.



## Blue plugged personal stereo V lead

Available in black or white in the following lengths:

	Black	White
600mm	DPV600	DPV600W
800mm	DPV800	DPV800W
1metre	DPV1M	DPV1MW
1.5metre	DPV1M5	DPV1M5W





## Blue plugged personal stereo single leads

Available in black or white in the following lengths:

800mm	DPS800	DPS800W
1metre	DPS1M	DPS1MW



The three pin plugs go into the direct input shoe appropriate for your hearing aids.



*Different hearing aids require different shoes – see Connevans Catalogue, section 2.*

*The blue plug will fit anywhere that a pair of headphones can be plugged in.*

**Also see Connevans Catalogue, section 11 – Listening for Pleasure**



**Connevans Limited**

Bridge House, 1 Nutfield Road, Merstham, Surrey RH1 3EB

Telephone: 01737 247571

info@connevans.com

Shop online: [www.DeafEquipment.co.uk](http://www.DeafEquipment.co.uk)